

1. Why is it necessary to correct the readings of the respirometers containing seeds with the readings taken from respirometers containing only glass beads? Your answer should refer to the concepts derived from the general gas law: $PV = nRT$ Where P = pressure of the gas V = volume of the gas n = number of moles of the gas R = the gas constant (its value is fixed) T = temperature of the gas.

It is necessary to correct readings based off the respirometer containing only glass beads because atmospheric temperature and pressure can change, even after the equilibration period. If the outside atmospheric pressure were to increase and everything else remained constant, $V = (nRT)/P$. This equation means that the measured volume of consumed O₂ would decrease when compared to the actual volume of consumed O₂ because the denominator P is becoming greater. If the temperature increased, $V = (nRT)/P$ holds true again. This equation means that the volume of consumed O₂ measured would be greater when compared to the actual volume of consumed O₂ because the numerator would be greater after multiplying the higher temperature by n and R. For both of these cases the opposite could also be the case, with a pressure decrease, an increase in measured volume of consumed O₂ when compared to the actual volume of consumed O₂ and with a temperature decrease a decrease in measured volume of consumed O₂ when compared to the actual volume of consumed O₂. To correct these changes in atmospheric conditions, a control respirometer will be used containing only glass beads. From this respirometer the other readings can be corrected.

2. What happens to the volume of the gas being measured (O₂ consumption or CO₂ production) when the temperature or pressure changes during the experiment? If pressure and temperature remain constant, will the volume of gas in the respirometers increase or decrease? Please explain. Hint: Several tutorials and animations explaining the general gas law are available online (e.g., <http://www.nclark.net/GasLaws>).

When atmospheric temperature and pressure change the O₂ consumed, what will be measured in this experiment, will change corresponding to $V = (nRT)/P$ as indicated in the table below:

atmospheric condition:	pressure increase	pressure decrease	temperature increase	temperature decrease
result on measure volume of O ₂ consumed when compared to actual volume of consumed O ₂ :	decrease	increase	increase	decrease

However if the atmospheric temperature and pressure remain constant, the volume of consumed O₂ will not change due to atmospheric conditions. The volume of consumed O₂ will still decrease as the experiment is being conducted as the germinating seeds use O₂. That consumption is what this experiment is designed to test for.